Finalidade:

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3. Baseman JB, Tully JG. *Mycoplasmas*: sophisticated, reemerging, and burdened by their notoriety. *Emerg Infect Dis* 1997;3:21.

4. Baseman JB, Reddy SP, Dallo SF. Interplay between mycoplasma surface proteins, airway cells, and the protean manifestations of mycoplasma-mediated human infections. *Am J Respir Crit Care Med* 1996;154(4 Pt 2): S137-44.

5. Berger RP, Wadowksy RM. Rhabdomyolysis associated with infection by *Mycoplasma pneumoniae*: a case report. *Pediatrics* 2000;105:433.

6. Bitnun A, Ford-Jones EL, Petric M, *et al.* Acute childhood encephalitis and *Mycoplasma pneumoniae*. *Clin Infect Dis* 2001;32:1674.

7. Bitnun A, Ford-Jones E, Blaser S, *et al.* *Mycoplasma pneumoniae* encephalitis. *Semin Pediatr Infect Dis* 2003;14: 96.

8. Blasco Patino F, Perez Maestu R, Lopez de Letona JM. Mechanisms of disease in *Mycoplasma pneumoniae* infection: clinical manifestations and complications. *Rev Clin Esp* 2004;204:365.

9. Broughton RA. Infections due to *Mycoplasma pneumoniae* in childhood. *Pediatr Infect Dis* 1986;5:71.

14. Cartner SC, Lindsey JR, Gibbs-Erwin J, *et al.* Roles of innate and adaptive immunity in respiratory mycoplasmosis. *Infect Immun* 1998;66:3485.

10. Cherry JD. Anemia and mucocutaneous lesions due to *Mycoplasma pneumoniae* infections. *Clin Infect Dis* 1993;17 Suppl 1:S47.

11. Chiba H, Pattanjitvilai S, Mitsuzawa H, *et al.* Pulmonary surfactant proteins A and D recognize lipid ligands on *Mycoplasma pneumoniae* and markedly augment the innate immune response to the organism. *Chest* 2003;123(3 Suppl):426S.

12. Christie LJ, Honarmand S, Talkington DF, *et al.* Pediatric encephalitis: what is the role of *Mycoplasma pneumoniae*? *Pediatrics* 2007;120:305-13.

13. Chryssanthopoulos C, Eboriadou M, Monti K, *et al.* Fatal disseminated intravascular coagulation caused by *Mycoplasma pneumoniae*. *Pediatr Infect Dis J* 2001;20:634.

14. Chua KB, Ngeow YF, Ng KB, *et al.* *Ureaplasma urealyticum* and *Mycoplasma hominis* isolation from cervical secretions of pregnant women and nasopharyngeal secretions of their babies at delivery. *Singapore Med J* 1998;39:300.

15. Clyde WA Jr. Clinical overview of typical *Mycoplasma pneumoniae* infections. *Clin Infect Dis* 1993;17 Suppl 1:S32.

16. Dallo SF, Baseman JB. Intracellular DNA replication and long-term survival of pathogenic mycoplasmas. *Microb Pathog* 2000;29:301.

17. Dominguez SR, Littlehorn C, Nyquist AC. *Mycoplasma hominis* endocarditis in a child with a complex congenital heart defect. *Pediatr Infect Dis J* 2006;25:851-2.

18. Donders GG, Van Bulck B, Caudron J, *et al.* Relationship of bacterial vaginosis and mycoplasmas to the risk of spontaneous abortion. *Am J Obstet Gynecol* 2000; 183:431.

19. Esposito S, Bosis S, Cavagna R, *et al.* Characteristics of *Streptococcus pneumoniae* and atypical bacterial infections in children 2-5 years of age with community-acquired pneumonia. *Clin Infect Dis* 2002;35:1345.

20. Esposito S, Blasi F, Droghetti R, *et al.* Double-blind, randomized, placebo controlled trial of clarithromycin therapy in children with recurrent acute wheezing [abstract G-1542]. In: 43rd Interscience Conference on Antimicrobial Agents and Chemotherapy; 2003 Sept 14-17; Chicago. Washington, DC: American Society for Microbiology; p. 288.

21. Ferwerda A, Moll HA, de Groot R. Respiratory tract infections by *Mycoplasma pneumoniae* in children: a review of diagnostic and therapeutic measures. *Eur J Pediatr* 2001;160:483.

22. File TM Jr, Tan JS, Plouffe JF. The role of atypical pathogens: *Mycoplasma pneumoniae*, *Chlamydia pneumoniae*, and *Legionella pneumophila* in respiratory infection. *Infect Dis Clin North Am* 1998;12:569-92, vii.

23. Fonseca-Aten M, Rios AM, Mejias A, *et al.* Treatment of experimental chronic pulmonary mycoplasmosis. *Int J Antimicrob Agents* 2006;28:253-8.

24. Fonseca-Aten M, Rios AM, Mejias A, *et al.* *Mycoplasma pneumoniae* induces host-dependent pulmonary inflammation and airway obstruction in mice. *Am J Respir Cell Mol Biol* 2005;32:201.

25. Fonseca-Aten M, Rios AM, Hardy RD. The role of atypical agents in respiratory tract disease. In: Kimpen JLL, Ramilo O, editors. The microbe-host interface in respiratory tract infections: current research and hot topics. Wymondham (UK): Horizon Bioscience; 2005. p. 133.
26. Foy HM. Infections caused by *Mycoplasma pneumoniae* and possible carrier state in different populations of patients. *Clin Infect Dis* 1993;17 Suppl 1:S37.
27. Gleason PP. The emerging role of atypical pathogens in community-acquired pneumonia. *Pharmacotherapy* 2002; 22(1 Pt 2):2S–11S; discussion 30S–32S.
28. Gupta SK, Sarosi GA. The role of atypical pathogens in community-acquired pneumonia. *Med Clin North Am* 2001;85:1349.
29. Halm EA, Teirstein AS. Clinical practice: management of community-acquired pneumonia. *N Engl J Med* 2002; 347:2039.
30. Hansbro PM, Beagley KW, Horvat JC, *et al.* Role of atypical bacterial infection of the lung in predisposition/ protection of asthma. *Pharmacol Ther* 2004;101:193.
31. Hardy RD, Atkinson TP, Cassell GH. Immune responses to *Mycoplasma*. In: Mestecky J, Lamm EL, Bienenstock J, *et al.* editors. *Mucosal immunology*. 3rd ed. New York: Elsevier; 2005. p. 1451–64.
32. Hardy RD, Ramilo O. *Mycoplasma* infections. In: Remington JS, Klein JO, Wilson CB, Baker CJ, editors. *Infectious diseases of the fetus and newborn infant*. 6th ed. Philadelphia, PA: Elsevier Saunders; 2006. p. 499–512.
33. Hardy RD, Jafri HS, Olsen K, *et al.* *Mycoplasma pneumoniae* induces chronic respiratory infection, airway hyperreactivity, and pulmonary inflammation: a murine model of infection-associated chronic reactive airway disease. *Infect Immun* 2002;70:649.
34. Himmelreich R, Plagens H, Hilbert H, *et al.* Comparative analysis of the genomes of the bacteria *Mycoplasma pneumoniae* and *Mycoplasma genitalium*. *Nucleic Acids Res* 1997;25:701.
35. Hirshberg SJ, Charles RS, Ettinger JB. Pediatric priapism associated with *Mycoplasma pneumoniae*. *Urology* 1996;47: 745.
36. Kaufmann A, Muhlradt PF, Gemsa D, *et al.* Induction of cytokines and chemokines in human monocytes by *Mycoplasma fermentans*-derived lipoprotein MALP-2. *Infect Immun* 1999;67:6303.
37. Kahane I. Oxidative damage induced by mycoplasmas. In: Razin S, Tully JG, editors. *Molecular and diagnostic procedures in mycoplasmaology*. New York: Academic Press; 1995. p. 415.
38. Kannan TR, Baseman JB. ADP-ribosylating and vacuolating cytotoxin of *Mycoplasma pneumoniae* represents unique virulence determinant among bacterial pathogens. *Proc Natl Acad Sci U S A* 2006;103:6724–9.
39. Kenney RT, Li JS, Clyde WA Jr, *et al.* *Mycoplasma* pericarditis: evidence of invasive disease. *Clin Infect Dis* 1993;17 Suppl 1:S58.
40. Kraft M, Cassell GH, Pak J, *et al.* *Mycoplasma pneumoniae* and *Chlamydia pneumoniae* in asthma: effect of clarithromycin. *Chest* 2002;121:1782.
41. Krause DC. *Mycoplasma pneumoniae* cytoadherence: unravelling the tie that binds. *Mol Microbiol* 1996;20:247.
42. Krause DC. *Mycoplasma pneumoniae* cytoadherence: organization and assembly of the attachment organelle. *Trends Microbiol* 1998;6:15.
43. Lieberman D, Schlaeffer F, Boldur I, *et al.* Multiple pathogens in adult patients admitted with community-acquired pneumonia: a one year prospective study of 346 consecutive patients. *Thorax* 1996;51:179.
44. Lo SC. *Mycoplasmas* and AIDS. In: Maniloff J, editor. *Mycoplasmas: molecular biology and pathogenesis*. Washington (DC): American Society for Microbiology; 1992. p. 525.
45. Luhrmann A, Deiters U, Skokowa J, *et al.* In vivo effects of a synthetic 2-kilodalton macrophage-activating lipopeptide of *Mycoplasma fermentans* after pulmonary application. *Infect Immun* 2002;70:3785.
46. Martin RJ, Kraft M, Chu HW, *et al.* A link between chronic asthma and chronic infection. *J Allergy Clin Immunol* 2001;107:595.
47. McCormack WM. Susceptibility of mycoplasmas to antimicrobial agents: clinical implications. *Clin Infect Dis* 1993;17 Suppl 1:S200.
48. Michelow IC, Olsen K, Lozano J, *et al.* Diagnostic utility and clinical significance of naso- and oropharyngeal samples used in a PCR assay to diagnose *Mycoplasma pneumoniae* infection in children with community-acquired pneumonia. *J Clin Microbiol* 2004;42:3339.
49. Miyashita N, Fukano H, Okimoto N, *et al.* Clinical presentation of community-acquired *Chlamydia pneumoniae* pneumonia in adults. *Chest* 2002;121:1776.
50. Miyashita N, Fukano H, Yoshida K, *et al.* Is it possible to distinguish between atypical pneumonia and bacterial pneumonia? Evaluation of the guidelines for community-acquired pneumonia in Japan. *Respir Med* 2004;98:952.
51. Miyashita N, Obase Y, Ouchi K, *et al.* Clinical features of severe *Mycoplasma pneumoniae* pneumonia in adults admitted to an intensive care unit. *J Med Microbiol* 2007;56:1625–9.
52. Mushegian AR, Koonin EV. A minimal gene set for cellular life derived by comparison of complete bacterial genomes. *Proc Natl Acad Sci U S A* 1996;93:10268.
53. Narita M, Yamada S. Two distinct patterns of central nervous system complications due to *Mycoplasma pneumoniae* infection. *Clin Infect Dis* 2001;33:916.
54. Niederman MS, Mandell LA, Anzueto A, *et al.* Guidelines for the management of adults with community-acquired pneumonia: diagnosis, assessment of severity, antimicrobial therapy, and prevention. *Am J Respir Crit Care Med* 2001;163:1730.
55. Nilsson AC, Bjorkman P, Persson K. Polymerase chain reaction is superior to serology for the diagnosis of acute *Mycoplasma pneumoniae* infection and reveals a high rate of persistent infection. *BMC Microbiol* 2008;8:93.
56. Orlicek SL, Walker MS, Kuhls TL. Severe *Mycoplasma pneumoniae* in young children with Down syndrome. *Clin Pediatr (Phila)* 1992;31:409.
57. Petitjean J, Vabret A, Gouarin S, *et al.* Evaluation of four commercial immunoglobulin G (IgG)- and IgM-specific enzyme immunoassays for diagnosis of *Mycoplasma pneumoniae* infections. *J Clin Microbiol* 2002;40:165.
58. Pfausler B, Engelhardt K, Kampf A, *et al.* Post-infectious central and peripheral nervous system diseases complicating *Mycoplasma pneumoniae* infection: report of three cases and review of the literature. *Eur J Neurol* 2002;9:93.
59. Potts JM, Ward AM, Rackley RR. Association of chronic urinary symptoms in women and *Ureaplasma urealyticum*. *Urology* 2000;55:486.
60. Rios AM, Mejias A, Chavez-Bueno S, *et al.* The impact of cethromycin (ABT-773) therapy on microbiologic, histologic, immunologic and respiratory indices in a murine model of *Mycoplasma pneumoniae* lower respiratory infection. *Antimicrob Agents Chemother* 2004;48:2897.
61. Salvatore CM, Fonseca-Aten M, Katz K, *et al.* Intranasal IL-12 therapy inhibits *Mycoplasma pneumoniae* clearance and sustains airway obstruction in murine pneumonia. *Infect Immun* 2008;76:732.
62. Seya T, Matsumoto M. A lipoprotein family from *Mycoplasma fermentans* confers host immune activation through Toll-like receptor 2. *Int J Biochem Cell Biol* 2002;34:901.
63. Sexually transmitted diseases treatment guidelines 2002. Centers for Disease Control and Prevention. *MMWR Recomm Rep* 2002;51(RR-6):1.
64. Sielaff TD, Everett JE, Shumway SJ, *et al.* *Mycoplasma hominis* infections occurring in cardiovascular surgical patients. *Ann Thorac Surg* 1996;61:99.
65. Smith R, Eviatar L. Neurologic manifestations of *Mycoplasma pneumoniae* infections: diverse spectrum of diseases: a report of six cases and review of the literature. *Clin Pediatr (Phila)* 2000;39:195.
66. Stacey CM, Munday PE, Taylor-Robinson D, *et al.* A longitudinal study of pelvic inflammatory disease. *Br J Obstet Gynaecol* 1992;99:994.
67. Stamm WE, Hicks CB, Martin DH, *et al.* Azithromycin for empirical treatment of the nongonococcal urethritis syndrome in men. *JAMA* 1995;274:545.
68. Stephan JL, Galambun C, Pozzetto B, *et al.* Aplastic anemia after *Mycoplasma pneumoniae* infection: a report of two cases. *J Pediatr Hematol Oncol* 1999;21:299.
69. Suzuki S, Yamazaki T, Narita M, *et al.* Clinical evaluation of macrolide-resistant *Mycoplasma pneumoniae*. *Anti microb Agents Chemother* 2006;50:709–12.

70. Szymanski M, Petric M, Saunders FE, *et al.* Mycoplasma pneumoniae pericarditis demonstrated by polymerase chain reaction and electron microscopy. Clin Infect Dis 2002;34:E16.

71. Tagliabue C, Salvatore CM, Techasaensiri C, *et al.* The impact of steroids given with macrolide therapy on experimental Mycoplasma pneumoniae respiratory infection. J Infect Dis 2008;198:1180–8.

72. Takeuchi O, Kaufmann A, Grote K, *et al.* Cutting edge: preferentially the R-stereoisomer of the mycoplasmal lipopeptide macrophage-activating lipopeptide-2 activates immune cells through a Toll-like receptor 2- and MyD88-dependent signaling pathway. J Immunol 2000; 164:554.

73. Tamura A, Matsubara K, Tanaka T, *et al.* Methylprednisolone pulse therapy for refractory Mycoplasma pneumoniae pneumonia in children. J Infect 2008;57:223–8.

74. Tarp B, Jensen JS, Ostergaard L, *et al.* Search for agents causing atypical pneumonia in HIV-positive patients by inhibitor-controlled PCR assays. Eur Respir J 1999;13:175.

75. Tay YK, Huff JC, Weston WL. Mycoplasma pneumoniae infection is associated with Stevens-Johnson syndrome, not erythema multiforme (von Hebra). J Am Acad Dermatol 1996;35:757.

76. Taylor-Robinson D. Mycoplasma genitalium: an update. Int J STD AIDS 2002;13:145.

77. Taylor-Robinson D, Davies HA, Sarathchandra P, *et al.* Intracellular location of mycoplasmas in cultured cells demonstrated by immunocytochemistry and electron microscopy. Int J Exp Pathol 1991;72:705.

78. Taylor-Robinson D. Infections due to species of Mycoplasma and Ureaplasma: an update. Clin Infect Dis 1996; 23:671.

79. Taylor-Robinson D, Gilroy CB, Thomas BJ, *et al.* Mycoplasma genitalium in chronic non-gonococcal urethritis. Int J STD AIDS 2004;15:21.

80. Templeton KE, Scheltinga SA, Graffelman AW, *et al.* Comparison and evaluation of real-time PCR, real-time nucleic acid sequence-based amplification, conventional PCR, and serology for diagnosis of Mycoplasma pneumoniae. J Clin Microbiol 2003;41:4366.

81. Totten PA, Schwartz MA, Sjostrom KE, *et al.* Association of Mycoplasma genitalium with nongonococcal urethritis in heterosexual men. J Infect Dis 2001;183:269.

82. Waites KB, Talkington DF. Mycoplasma pneumoniae and its role as a human pathogen. Clin Microbiol Rev 2004;17: 697.

83. Waites KB, Crabb DM, Bing X, *et al.* In vitro susceptibilities to and bactericidal activities of garenoxacin (BMS284756) and other antimicrobial agents against human mycoplasmas and ureaplasmas. Antimicrob Agents Chemother 2003;47:161.

84. Waites KB, Crabb DM, Duffy LB. In vitro activities of ABT773 and other antimicrobials against human mycoplasmas. Antimicrob Agents Chemother 2003;47:39.

85. Winchell JM, Thurman KA, Mitchell SL, *et al.* Evaluation of three real-time PCR assays for detection of Mycoplasma pneumoniae in an outbreak investigation. J Clin Microbiol 2008;46:3116–8.

86. Winner F, Rosengarten R, Citti C. In vitro cell invasion of Mycoplasma gallisepticum. Infect Immun 2000;68:4238.

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