Finalidade:

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ATCC®: American Type Culture Collection

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1. Baker, C.N., C. Thornsberry, and R.W. Hawkinson. 1983. Inoculum standardization in antimicrobial susceptibility testing: evaluation of overnight agar cultures and the rapid inoculum standardization system. *J. Clin. Microbiol.* 17:450-457.
2. Barry, A.L., F. Garcia, and L.D. Thrupp. 1970. An improved single-disk method for testing the antibiotic susceptibility of rapidly-growing pathogens. *Am. J. Clin. Pathol.* 53:149-158.
3. Bauer, A.W., W.M.M. Kirby, J.C. Sherris, and M. Turck. 1966. Antibiotic susceptibility testing by a standardized single disk method. *Am. J. Clin. Pathol.* 45:493-496.
4. CLSI. Approved standard: M02-A12 - Performance standards for antimicrobial disk susceptibility tests. CLSI, Wayne, PA, USA. Search for latest version at www.clsi.org.
5. CLSI - Performance Standards for Antimicrobial Susceptibility. CLSI supplement M100S. Wayne, PA: Clinical and Laboratory Standards Institute. Search for latest version at www.clsi.org.
6. CLSI. Performance Standards for Antimicrobial Susceptibility Testing. 26th ed. CLSI supplement M100S. Wayne, PA: Clinical and Laboratory Standards Institute; 2016.
7. CLSI. Approved standard: M7. Methods for dilution antimicrobial susceptibility tests for bacteria that grow aerobically. CLSI, Wayne, PA, USA. Search for latest version at www.clsi.org.
8. Data on file. Becton Dickinson.
9. D'Amato, R.F., and C. Thornsberry. 1979. Calcium and magnesium in Mueller-Hinton agar and their influence on disk diffusion susceptibility results. *Current Microbiol.* 2:135-138.
10. Ericsson, H.M., and J.C. Sherris. 1971. Antibiotic sensitivity testing. Report of an international collaborative study. *Acta Pathol. Microbiol. Scand. Sec. B, Suppl.* 217.
11. EUCAST Disk Diffusion Method for Antimicrobial Susceptibility Testing. Search for latest version at <http://www.eucast.org>.
12. Ferone, R., S.R.M. Bushby, J.J. Burchall, W.D. Moore, and D. Smith. 1975. Identification of Harper-Cawston factor as thymidine phosphorylase and removal from media of substances interfering with susceptibility testing to sulfonamides and diaminopyrimidines. *Antimicrob. Agents Chemother.* 7:91-98.
13. Haltiner, R.C., P.C. Migneault, and R.G. Robertson. 1980. Incidence of thymidine-dependent enterococci detected on Mueller-Hinton agar with low thymidine content. *Antimicrob. Agents Chemother.* 18:365-368.
14. Hindler, J.A., and C.B. Anderbied. 1985. Effect of the source of Mueller-Hinton agar and resistance frequency on the detection of methicillin-resistant *Staphylococcus aureus*. *J. Clin. Microbiol.* 21:205-210.
15. Hollick, G.E., and J.A. Washington II. 1976. Comparison of direct and standardized disk diffusion susceptibility testing of urine cultures. *Antimicrob. Agents Chemother.* 9:804-809.
16. Jorgensen, J.H., and J.D. Turnidge. 2003. Susceptibility test methods: dilution and disk diffusion methods. In: Murray, P. R., E. J.

Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.

17. Johnson, J.E., and J.A. Washington II. 1976. Comparison of direct and standardized antimicrobial susceptibility testing of positive blood cultures. *Antimicrob. Agents Chemother.* 10:211-214.

18. Koch, A.E., and J.J. Burchall. 1971. Reversal of the antimicrobial activity of trimethoprim by thymidine in commercially prepared media. *Appl. Microbiol.* 22:812-817.

19. Matuschek, E., D.F. Brown, and G. Kahlmeter. Development of the EUCAST disk diffusion antimicrobial susceptibility testing method and its implementation in routine microbiology laboratories. *Clin. Microbiol. Infect.* 2014; 20 (4): 255-66.

20. Maskell, R., O.A. Okubadejo, R.H. Payne, and L. Peard. 1977. Human infections with thymine-requiring bacteria. *J. Med. Microbiol.* 11:33-45.

21. Murray, B.E. 1990. The life and times of the Enterococcus. *Clin. Microbiol. Rev.* 3:46-65.

22. Mueller, J.H., and J. Hinton. 1941. A protein-free medium for primary isolation of the gonococcus and meningococcus. *Proc. Soc. Exp. Biol. Med.* 48:330-333.

23. Neumann, M.A., D.F. Sahn, C. Thornsberry, J.E. McGowan, Jr. Cumitech 6A, New developments in antimicrobial agent susceptibility testing: a practical guide. Coordinating ed., J.E. McGowan, Jr. American Society of Microbiology, Washington, D.C. 1991.

24. Pollock, H.M., B.H. Minshew, M.A. Kenny, and F.D. Schoenknecht. 1978. Effect of different lots of Mueller-Hinton Agar on the interpretation of the gentamicin susceptibility of *Pseudomonas aeruginosa*. *Antimicrob. Agents Chemother.* 14:360-367.

25. Reller, L.G., F.D. Schoenknecht, M.A. Kenny, and J.C. Sherris. 1974. Antibiotic susceptibility testing of *Pseudomonas aeruginosa*: selection of a control strain and criteria for magnesium and calcium content in media. *J. Infect. Dis.* 130:454-463.

26. Ryan, K.J., F.D. Schoenknecht, and W.M.M. Kirby. 1970. Disc sensitivity testing. *Hospital Practice* 5:91-100.

27. Skov, R., A.R. Larsen, A. Kearns, M. Holmes, C. Teale, G. Edwards, and R. Hill. Phenotypic detection of *mecC*-MRSA: ceftiofuran is more reliable than oxacillin. 2014. *J. Antimicrob. Chemother.* 69:133-135.

28. The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. Search for latest version at <http://www.eucast.org>.

29. Thornsberry, C., T.L. Gavan, and E.H. Gerlach. 1977. Cumitech 6, New developments in antimicrobial agent susceptibility testing. Coordinating ed., J.C. Sherris. American Society of Microbiology, Washington, DC.

30. The European Committee on Antimicrobial Susceptibility Testing. Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST. Search for latest version at <http://www.eucast.org>.

31. The European Committee on Antimicrobial Susceptibility Testing. Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST. Version 6.0, 2016. <http://www.eucast.org>.

32. Washington, J.A., and G.L. Woods. 1995. Antimicrobial susceptibility tests: dilution and disk diffusion methods. P. 1327-1341. In Muarry, P.R., Baron, E.J., Pfaller, M.A., Tenover, F.C., and Tenover, R.H. (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C. 1995.

33. Waterworth, P.M., and M. Del Piano. 1976. Dependability of sensitivity tests in primary culture. *J. Clin. Pathol.* 29:179-184.

34. Wegner, D.L., C.R. Mathis, and T.R. Neblett. 1976. Direct method to determine the antibiotic susceptibility of rapidly growing blood pathogens. *Antimicrob. Agents Chemother.* 9:861-862.

35. Woods, G.L., and J.A. Washington. 1995. Antibacterial susceptibility tests: dilution and disk diffusion methods, p. 1327-1341. In P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.C. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, DC.

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