Finalidade:

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11. REFERÊNCIAS

1. Ajello, L., L.K. Georg, W. Kaplan, and L. Kaufman. 1963. CDC laboratory manual for medical mycology. PHS Publication No. 994, U.S. Government Printing Office, Washington, D.C.
2. Almeida-Paes, R. et al. Immunoglobulins G, M, and A against *Sporothrix schenckii* exoantigens in patients with sporotrichosis before and during treatment with itraconazole. *Clinical and Vaccine Immunology*. v. 14, n. 9, p. 1149-1157, 2007.
3. Attili, S. D. Importância e sistemática de fungos filamentosos. Campinas: Fundação Tropical de PesquisaTecnologia André Tosello, 1990.
4. Biswas S, Dijck PV, Datta A. Environmental sensing and signal transduction pathways regulating morphopathogenic determinants of *Candida albicans*. *Microbiol Mol Biol R* 2007; 71(2): 348-376.
5. Brown AJB, Gow NAR. Regulatory networks controlling *Candida albicans* morphogenesis. *Trends Microbiol* 1999; 7(8): 333-338.
6. Crespo-Erchiga, V. ; Gmez-Moyano, E. ; C, M. La pitiriasis versicolor y las levaduras del gÉnero *Malassezia*. *Actas Dermo-Sifilograficas*. v. 99, n. 10, p. 764-771, 2008.
7. Calderone RA, Fonzi WA. Virulence factors of *Candida albicans*. *Trends Microbiol* 2001; 9(7): 327-335
8. Chaves GM, Cavalvanti MAQ, Porto ALF. Pathogenicity characteristics of stocked and fresh yeast strains. *Braz J Microbiol* 2003; 34: 197-202.
9. ChaffinWL, López-Ribot JL, Casanova M, Gozalbo D, Martínez JP. Cell wall and secreted proteins of *Candida albicans*: identification, function, and expression. *Microbiol Mol Biol R* 1998; 62(1): 130-180.
10. Colombo AL, Gimarões T. Epidemiologia das infecções hematogênicas por *Candida* spp. *Rev Soc Bras Med Trop* 2003; 36(5): 599-607.
11. Crampin H, Finley K, Geramid-Nejad M, Court H, Gale C, Berman J et al. *Candida albicans* hyphae have a Spitzenkörper that

is distinct from the polarisome found in yeast and pseudohyphae. J Cell Sci 2005; 118(13): 2935-2947.

12. Fromtling, R.A. 1995. Mycology. In: P.R. Murray, E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.

13. Forche A, May G, Beckerman J, Kauffman S, Becker J, Magee PT. A system for studying genetic changes in *Candida albicans* during infection. Fungal Genet Biol 2003; 39: 38-50.

14. Knight SAB, Vilaire G, Lesuisse E, Dancis A. Iron acquisition from transferrin by *Candida albicans* depends on the reductive pathway. Infect Immun 2005; 73(9): 5482-5492.

15. Haley, L.D., J. Trandel, and M.B. Coyle. 1980. Cumitech 11, Practical methods for culture and identification of fungi in the clinical microbiology laboratory. Coordinating ed., J.C. Sherris. American Society for Microbiology, Washington, D.C.

16. Hube B. From commensal to pathogen: stage- and tissue-specific gene expression of *Candida albicans*. Curr Opin Microbiol 2004; 7: 336-341.

17. Hromatka BS, Noble SM, Johnson AD. Transcriptional response of *Candida albicans* to nitric oxide and the role of the YHB1 gene in nitrosative stress and virulence. Mol Biol Cell 2005; 16: 4814-4826.

18. Jong AY, Stins MF, Huang S-H, Chen SHM, Kim KS. Traversal of *Candida albicans* across human blood-brain barrier in vitro. Infect Immun 2001; 69(7): 4536-4544.

19. Kwon-Chung, K.J., and J.E. Bennett. 1992. Medical mycology. Lea & Febiger, Philadelphia.

20. Larone, D.H. 2002. Medically important fungi: a guide to identification. 4th ed. American Society for Microbiology, Washington, D.C.

21. Lopes-Bezerra, L. M. ; Schubach, A. O. ; Costa, R. O. *Sporothrix schenckii* and sporotrichosis. Anais da Academia Brasileira de Ciências. v. 78, n. 2, p. 293-308, 2006.

22. MacFaddin, J.F. 1985. Media for isolation-cultivation-identification-maintenance of medical bacteria. vol. I. Williams & Wilkins, Baltimore.

23. Menezes EA, Guerra ACP, Rodrigues RCB, Peixoto MMLV, Lima LS, Cunha FA. Isolamento de *Candida* spp. no mamilo de lactantes do banco de leite humano da universidade federal do ceará e teste de suscetibilidade a antifúngicos. J Bras Patol Med Lab 2004; 40(5): 299-305.

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27. Odds FC. Effects of temperature on anti-*Candida* activities of antifungal antibiotics. Antimicrob Agents Ch 1993; 37(4): 685-691.

28. Pfaller MA, Diekema DA. Epidemiology of invasive candidiasis: a persistent public health problem. Clin Microbiol Rev 2007; 20(1): 133-163.

29. Richardson, M.D.; Warnock, D.W. - Fungal infection: diagnosis and management. Oxford: Blackwell, 1993

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32. Soll DR, Lockhart SR, Zhao R. Relationship between Switching and Mating in *Candida albicans*. Eukaryotic Cell 2003; 2(3): 390-397.

33. Sudbery P, Gow N, Berman J. The distinct morphogenic states of *Candida albicans*. Trends Microbiol 2004; 12(7): 313-324.

34. Sutton, D.A. 2003. Specimen collection, transport, and processing: mycology. In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.

35. Summerbell, R.C. 2003. Trichophyton, Microsporum, Epidermophyton, and agents of superficial mycoses. In: Murray, P. R., E. J. Baron, J.H. Jorgensen, M. A. Pfaller, and R. H. Tenover (ed.). Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.

36. Romani L, Bistoni F, Puccetti P. Adaptation of *Candida albicans* to the host environment: the role of morphogenesis in virulence and survival in mammalian hosts. Curr Opin Microbiol 2003; 6: 338-343.

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