## Finalidade:

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## 11. REFERÊNCIAS

1. Baker, C.N.; Thomsberry, C.; and Hawkinson, R.W. 1983. Inoculum standardization in antimicrobial susceptibility testing: evaluation of overnight agar cultures and the rapid inoculum standardization system. *J. Clin. Microbiol.* 17:450-457.
2. Barry, A.L.; Garcia, F.; and Thrupp, L.D. 1970. An improved single-disk method for testing the antibiotic susceptibility of rapidly-growing pathogens. *Am. J. Clin. Pathol.* 53:149-158.

3. CLSI. Approved standard: M02-A12 - Performance standards for antimicrobial disk susceptibility tests. CLSI, Wayne, PA, USA. Search for latest version at [www.clsi.org](http://www.clsi.org).
4. CLSI - Performance Standards for Antimicrobial Susceptibility. CLSI supplement M100S. Wayne, PA: Clinical and Laboratory Standards Institute. Search for latest version at [www.clsi.org](http://www.clsi.org).
5. CLSI - Approved standard: M7. Methods for dilution antimicrobial susceptibility tests for bacteria that grow aerobically. CLSI, Wayne, PA, USA. Search for latest version at [www.clsi.org](http://www.clsi.org).
6. D'Amato, R.F. and Thornsberry, L.D. 1979. Calcium and magnesium in Mueller-Hinton agar and their influence on disk diffusion susceptibility results. *Current Microbiol.* 2:135-138.
7. Disk Diffusion Method for Antimicrobial Susceptibility Testing. Search for latest version at <http://www.eucast.org>.
8. ISO: Susceptibility testing of infectious agents and evaluation of performance of antimicrobial susceptibility test-devices – Part 1: Broth micro-dilution reference method for testing the in vitro activity of antimicrobial agents against rapidly growing aerobic bacteria involved in infectious diseases. ISO 20776-1:2019(E).
9. Johnson, J.E. and Washington J.A II. 1976. Comparison of direct and standardized antimicrobial susceptibility testing of positive blood cultures. *Antimicrob. Agents Chemother.* 10:211-214.
10. Kavanagh, A.; Ramu, S.; Gong, Y.; Gong, M.A. and Blaskovich, M.A.T. Effects of microplate type and broth additives on microdilution MIC susceptibility assays. *Antimicrob Agents Chemother* 63:e01760-2018.
11. Koneman, E. W. Diagnóstico microbiológico. Guanabara Koogan 6ed., 2010.
12. Matuschek, E., D.F. Brown, and G. Kahlmeter. Development of the EUCAST disk diffusion antimicrobial susceptibility testing method and its implementation in routine microbiology laboratories. *Clin. Microbiol. Infect.* 2014; 20(4): 255-66.
13. Murray, B.E. 1990. The life and times of the *Enterococcus*. *Clin. Microbiol. Rev.* 3:46-65.
14. Mueller, J.H., and Hinton J. 1941. A protein-free medium for primary isolation of the *Gonococcus* and *Meningococcus*. *Proc. Soc. Exp. Biol. Med.* 48:330-333.
15. Oplustil, C.P.; Zoccoli, C.M.; Tobouti, N.R. e Scheffer, M.C. Procedimentos Básicos em Microbiologia Clínica. 4.ed. Sarvier, São Paulo, 2020.
16. The European Committee on Antimicrobial Susceptibility Testing (EUCAST). Breakpoint tables for interpretation of MICs and zone diameters. Search for latest version at <http://www.eucast.org>.
17. The European Committee on Antimicrobial Susceptibility Testing (EUCAST). Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST. Search for latest version at <http://www.eucast.org>.
18. Vaudaux, P.; Huggler, E.; Bernard, L.; Ferry, T.; Renzoni, A., and Lew, D.P. 2010. Underestimation of vancomycin and teicoplanin MICs by broth microdilution leads to underdetection of glycopeptide-intermediate isolates of *Staphylococcus aureus*. *Antimicrobial Agents and Chemotherapy*, 54(9):3861-3870.
19. Washington, J.A. and Woods, G.L. 1995. Antimicrobial susceptibility tests: dilution and disk diffusion methods. P. 1327-1341. In Muarry, P.R., Baron, E.J., Pfaller, M.A., Tenover, F.C., and Tenover, R.H. (ed.), *Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, D.C. 1995.
20. Waterworth, P.M. and Del Piano, M. 1976. Dependability of sensitivity tests in primary culture. *J. Clin. Pathol.* 29:179-184.

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